


Accessories

Intended to enable connections during the transition to ENFit® ISO80369-3

Tethered Cap

	30 pcs / box
Product Code:	ISO2UAD
NHSSC:	FWM4228



ISO2UAD Female universal funnel connector to female ENFit – will fit male universal connector on legacy products and convert to female ENFit

Transition Adaptors

- ENFit® compatible
- Latex free



	100 pcs / box
Product Code:	UAD2ISO
NHSSC:	FWM4223

UAD2ISO Male universal funnel connector to Male ENFit – will fit female universal connector on legacy products and convert to male ENFit



	10 pcs / box
Product Code:	TRN202
NHSSC:	FWM2436

Glow Green Connector

Transition Adaptors



Accessories

Double male ENFit® connector (ISO 80369-3), converts any ENFit® single port feeding tube to Y-Port, sterile

Y-Port Adaptor

- Converts any ENFit® single port feeding tube to Y-port
- Double male ENFit® connectors with caps
- ENFit® compliant



20 pcs / box

Product Code:	YPA-ISO
NHSSC:	FWM2432

Colour Coded Purple for Enteral Feeding

Winged Luer Adaptor

A 10 ml syringe with a Winged Luer Adapter is include with each AMT G-JET® Button and micro G-JET®. The Winged Luer Adapter facilitates direct G & J port access for physicians to inject contrast dye during placement/ replacement of the G-JET® device.

10 pcs / box

Product Code:	4-7000
NHSSC:	TBC

Adaptor

