

Enteral Syringe & Push Fit Cap

Enteral ISOSAF single use Enteral Syringe with female ENFit® connection (ISO 80369-3)

Low dose connectors are available

ENFit® Syringe & Push Fit Cap

- Sterile
- DEHP & latex free
- Clear barrel markings
- 60ml syringes have eccentric tips
- Rubber gasket
- Full thread
- Low dose option available in 1ml, 3ml and 6ml syringes
- Single use



Includes Push Fit Cap

1ml 100 pcs / box
Product Code: ELDISO-PFC-1




3ml 100 pcs / box
Product Code: ELDISO-PFC-3



6ml 100 pcs / box
Product Code: ELDISO-PFC-6




12ml 100 pcs / box
Product Code: EISO-PFC-12




20ml 100 pcs / box
Product Code: EISO-PFC-20



35ml 50 pcs / box
Product Code: EISO-PFC-35



60ml 50 pcs / box
Product Code: EISO-PFC-60



Enteral Syringe Single Use