


Accessories

Intended to extend Enteral feeding tubes or for use with Syringe drivers. Male & female ENFit® connectors (ISO 80369-3)



Extension Set

Size:	25 pcs / box Product Code:	NHSSC:
75cm	ESCF.075-ISOSAF	FWM2634
150cm	ESCF.150-ISOSAF	FWM2635

Enteral Extension Sets

- Bifurcated versions have priming volumes of less than 2ml
- ENFit® compliant for patient safety
- Slide clamp for one-handed use
- Latex free & DEHP free
- Choice of lengths



Bifurcated Extension Set

Size:	25 pcs / box Product Code:	NHSSC:
75cm	BESCF.075.ISOSAF	FWM85065
143cm	BESCF.143.ISOSAF	FWM85063

Enteral Extension Sets

