

Nutricare Flo

Nutricare Flo is an open-ended polyurethane tube with clearly numbered cm markings

ENFit® (ISO 80369-3) connector

Enteral Nutricare Flo

- Male ENFit® connector
- **DEHP & latex free**
- Made of PU
- **Radiopaque stripe**
- Clear cm markings
- **Includes fixation plaster**
- 2 non-opposing outlet holes
- **Rounded open distal end**
- *“Open ended feeding tubes have been reported to lead to a reduction in bacterial colonisation”¹*

¹ Enteral Feeding Tube Design and Differential Bacterial Overgrowth: An In Vitro Comparison. Amy Presti, MD, MSCR; Ruth Snyder, NNP, BC: Neonatal INTENSIVE CARE Vol. 26 No. 4 n July-August 2013



4Fr

25 pcs / box
Product Code:

50cm

NCF04.50

75cm

NCF04.75

6Fr

25 pcs / box
Product Code:

50cm

NCF06.50

60cm

NCF06.60

75cm

NCF06.75

10Fr

25 pcs / box
Product Code:

75cm

NCF10.75

5Fr

25 pcs / box
Product Code:

50cm

NCF05.50

60cm

NCF05.60

75cm

NCF05.75

8Fr

25 pcs / box
Product Code:

50cm

NCF08.50

60cm

NCF08.60

75cm

NCF08.75

Enteral Nutricare Flo