

Accessories

For keeping ENFit® connectors
& caps free of debris

Enteral Connector Brush



- ENFit® connectors have larger and deeper grooves, meaning increased area for debris to collect and potentially cause connection issues
- Using a brush after every use helps to prevent the accumulation of debris, and the feeding tube cap from adhering to the connector
- Keeping the connector and cap free of debris can help to extend the lifespan of the tube

In line with GEDSA ENFit cleaning Procedures 2018 recommendations
Ref <http://stayconnected.org/wpcontent/uploads/2018/05/ENFit-Cleaning-Procedures-All-Tubes.pdf>

50 pcs / box

Product Code: **ECB-NS**

NHSSC: **FWM3176**

